

Maintenance Log Sheet – ABL90FLEX

HOSPITAL: _____ LOCATION: _____ MONTH/YEAR _____

Lab/Ward staff must fill and initial the maintenance log each day

Daily	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Check AQUIRE connectivity																															
Check Analyser Status-if traffic light is not green, take appropriate action as per analyser need																															
Replace consumables as per analyser needs (Solution pack, sensor cassette, paper roll)																															
Check calibration and QC status, troubleshoot as per analyser needs																															
**Print and enter QC's in Pathnet (site specific)																															
Monthly																															
Run QAP																															
Quarterly																															
Perform tHb Calibration																															
Replace Inlet probe & Connection gasket																															
Initial																															

Issue Log

Date	Operator	Description of Problem/ Error Code	Action Taken	C.A.R Raised (if yes, C.A.R #)	Initial